



## PRODUCT SPECIFICATION

|                      |  |
|----------------------|--|
| catalog number:      | <b>A66013</b>  |
| product description: | ACTH Monoclonal Antibody, Purified   |
| lot number:          | [lot specific]   |
| immunogen:           | ACTH conjugate BSA   |
| source:              | Murine   |
| form:                | Liquid   |
| purification method: | Protein A chromatography   |
| purity:              | >98%   |
| buffer:              | 50mM NaPO <sub>4</sub> , 0.9% NaCl, pH 7.3   |
| isotype:             | IgG <sub>1</sub>   |
| preservative:        | 0.1% Proclin 300   |
| concentration:       | [lot specific] mg/mL   |
| specificity:         | This antibody is specific for human CLIP. It does not react with the N-terminal portion of ACTH. This is the detection antibody and can be used for conjugation and pairs with A66010. |
| application/notes:   | RIA, EIA, immunohistochemistry, and western blot   |
| storage:             | 2-8°C  |
| hazard/biohazard:    | There is no known biohazard associated with this product.  |