



PRODUCT SPECIFICATION

catalog number:	A66013
product description:	ACTH Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	ACTH conjugate BSA
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>98%
buffer:	PBS, pH 7.2±0.2, 0.02% Proclin
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
cross-reactivity:	n/a
application/notes:	RIA, EIA, immunohistochemistry, and western blot
storage:	2-8°C
hazard/biohazard:	There is no known biohazard associated with this product.