



PRODUCT SPECIFICATION

catalog number:	A66014
product description:	ACTH Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human ACTH (18-39)
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1 M BBS, pH 8
isotype:	IgG ₁ , Kappa
concentration:	[lot specific] mg/mL
specificity:	Recognizes Human ACTH (18-39) and Human ACTH (1-39). No cross reactivity to ACTH(1-24)
expiration date:	[lot specific]
storage:	2-8°C