



## PRODUCT SPECIFICATION

catalog number:	<b>A66022</b>
product description:	Osteocalcin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human recombinant Osteocalcin
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1M PBS, pH 7.4
isotype:	IgG <sub>1</sub> , Kappa
concentration:	[lot specific] mg/mL
specificity:	Recognizes native and recombinant Human Osteocalcin
pairing:	This antibody can be used as the detection antibody in ELISA.  Pair with A66021
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	There is no known biohazard associated with this product.