



PRODUCT SPECIFICATION

CATALOG NUMBER: A66030

PRODUCT:	Human Parathyroid Hormone (PTH N-terminal) Monoclonal Antibody
SOURCE:	Tissue Culture Supernatant
IMMUNOGEN:	Synthetic PTH (aa1-34)
PURIFICATION:	Protein A affinity chromatography
BUFFER:	Phosphate buffered saline, pH 7.2
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8 degrees short term, -20 degrees long term
SPECIFICITY:	This antibody reacts with 1-34 PTH
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.