



PRODUCT SPECIFICATION

catalog number:	A66051
product description:	PTH (C-Terminal) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant fragment of C-terminal human PTH protein
source:	Murine; in vitro production
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN ₃
isotype:	IgG _{2b} , kappa
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for human PTH.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.