



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A68040**

PRODUCT:	Carbamazepine Monoclonal Antibody
HOST:	Mouse
IMMUNOGEN:	Carbamazepine-BSA
PURITY:	Protein G chromatography
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	None
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8°C short term (one month), -20°C long term, Avoid repeated freeze thaw cycles
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.