



PRODUCT SPECIFICATION

catalog number:	A71040
product description:	Adiponectin/Acrp30 Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant human Adiponecin (Glu19-Asn244)
source:	Murine; cell culture supernatant
form:	Liquid
purification method:	Protein A/G chromatography; 0.2µm filtered
buffer:	PBS, pH 7.4
preservative:	None
isotype:	IgG _{2b}
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody Detects an epitope in the globular domain (aa 104-244) of human Adiponectin.
cross-reactivity:	This antibody shows no cross-reactivity with recombinant mouse Adiponectin or rh6Ckine.
storage:	Upon receipt store -20°C (up to 12 months); After opening: 2-8°C (one month); -20°C (up to 6 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.