



PRODUCT SPECIFICATION

catalog number:	A71042
product description:	Adiponectin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant human Adiponectin HEK293
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
preservative:	None
isotype:	IgG _{2b, kappa}
concentration:	[lot specific] mg/mL
specificity:	recognizes recombinant Human Adiponectin
storage:	2-8°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.