



PRODUCT SPECIFICATION

catalog number:	A72061
product description:	Coronavirus Nucleoprotein (NP)
lot number:	[lot specific]
source:	Murine ascites
form:	liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Reactive with the NP of native SARS and MERS coronaviruses. Non reactive with Coronavirus strains 22E, NL63 and OC43
notes/applications:	ELISA, IFA and Western Blot.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.