



## PRODUCT SPECIFICATION

catalog number:	<b>A72063</b>
product description:	Coronavirus Nucleoprotein (NP)
lot number:	[lot specific]
source:	Murine ascites
form:	liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	Reactive with the NP of native SARS and MERS coronaviruses.  Non reactive with Coronavirus strains 22E, NL63 and OC43
notes/applications:	ELISA, IFA and Western Blot.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.