



PRODUCT SPECIFICATION

CATALOG NUMBER: A78040

PRODUCT:	Beta-2-Microglobulin Monoclonal Antibody
IMMUNOGEN:	Human Beta-2-Microglobulin
FORM:	Purified IgG from mouse ascites
PURIFICATION:	Protein A chromatography
PURITY:	98% by SDS-PAGE
BUFFER:	50 mM Sodium Phosphate, 500 mM Sodium Chloride, 10 % glycerol (w/v), pH 8.0
PRESERVATIVE:	0.05% sodium azide
SUBCLASS:	IgG1, Kappa
STORAGE CONDITIONS:	2-8°C, Do not freeze
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.