



PRODUCT SPECIFICATION

catalog number:	A78052
product description:	β -2 Microglobulin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human β -2 microglobulin
source:	Murine ascites
form:	Liquid
purification method:	Affinity chromatography
buffer:	15mM phosphate, 0.15M NaCl, pH 7.2
preservative:	0.09% NaN ₃
isotype:	IgG _{1, kappa}
concentration:	[lot specific] mg/mL
specificity:	Human β -2 microglobulin
pairing:	This antibody can be used as the detection antibody in a sandwich EIA and paired with A78050 as the capture antibody.
storage:	2-8°C short term; -70°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.