



PRODUCT SPECIFICATION

catalog number:	A82260
product description:	Cyfra21-1 Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	MCF-7 cell
source:	Murine ascites
form:	Liquid
purity:	>90%
buffer:	10mM PBS, 150mM NaCl, pH 7.4
preservative:	0.1% P300
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL
pairing:	This antibody can be used in sandwich ELISA paired with either Cyfra21-1 MAb A82261 or HRP-Cyfra21-1 MAb T82261.
expiration date:	[lot specific]
storage:	2-8°C (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.