



PRODUCT SPECIFICATION

catalog number:	A82262
product description:	Cyfra21-1 Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Cell Culture
form:	Liquid
purity:	>90%
buffer:	10mM PBS, 150mM NaCl, pH 7.4
preservative:	0.1% P300
concentration:	[lot specific] mg/mL
pairing:	This antibody can be used in as capture in sandwich ELISA paired with either A82263 or HRP-Cyfra21-1 MAb T82261.
expiration date:	[lot specific]
storage:	2-8°C (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.