



## PRODUCT SPECIFICATION

|                      |  |
|----------------------|--|
| catalog number:      | <b>A82263</b>  |
| product description: | Cyfra21-1 Recombinant Antibody   |
| lot number:          | [lot specific]   |
| source:              | CHO Cells  |
| form:                | Liquid   |
| purity:              | >90%   |
| buffer:              | 10mM PBS, 150mM NaCl, pH 7.4   |
| preservative:        | 0.1% P300  |
| concentration:       | [lot specific] mg/mL   |
| pairing:             | This antibody can be used in as conjugate in sandwich ELISA paired with A82262 |
| expiration date:     | [lot specific]   |
| storage:             | 2-8°C (two weeks); -20°C long term avoiding multiple freeze thaw cycles.       |
| hazard/biohazard:    | There is no known hazard or biohazard associated with this product.            |