



PRODUCT SPECIFICATION

catalog number:	A83079601P
product description:	Chromogranin A Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1M PBS, pH 7.4
preservative:	0.5g/L NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.