



PRODUCT SPECIFICATION

CATALOG NUMBER: A83200

PRODUCT:	Intrinsic Factor Monoclonal Antibody
IMMUNOGEN:	Intrinsic Factor from porcine gastric mucosa
PURITY RESULTS:	HPLC/FPLC IgG content > 95%
BUFFER:	0.05M sodium acetate, 0.05M TRIS, 0.5 M sodium chloride, pH 7.0
PRESERVATIVE:	0.02% sodium azide
APPEARANCE:	Clear solution; final product .22u filtered
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
PRODUCTION:	This antibody was produced in vitro in serum-free media.
APPLICATION:	EIA
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.