



PRODUCT SPECIFICATION

catalog number:	A86154
product description:	CD45 Leucocyte Marker Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Stimulated human leukocytes
source:	Murine; in vitro production
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN ₃
isotype:	IgG _{2b} , kappa
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for human CD45 and is known to work in flow cytometry.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.