



PRODUCT SPECIFICATION

catalog number:	A86156
product description:	CD45 Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human peripheral blood mononuclear cells
source:	Murine; <i>in vitro</i> production
form:	Liquid
purification method:	Protein A or G chromatography; filtered in 0.2 µm solution
buffer:	PBS, pH 7.4
preservative:	none
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for human CD45 and recognizes all isoforms of human CD45.
storage:	Upon receipt store -20°C (up to 12 months); After opening: 2-8°C (one month); -20°C (up to 6 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.