



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A87001535P |
| product description: | NFL-Ct (Neurofilament Light Chain) Monoclonal Antibody |
| lot number: | [lot specific] |
| source: | Murine |
| form: | Liquid |
| purity: | >95% (SDS PAGE) |
| buffer: | Phosphate Buffered Saline (PBS), pH 7.4 |
| preservative: | 0.1% NaN ₃ |
| isotype: | IgG ₁ |
| concentration: | [lot specific] mg/mL |
| expiration date: | [lot specific] |
| storage: | 2-8°C |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |