



PRODUCT SPECIFICATION

catalog number:	A87002535P
product description:	NFL-Ct (Neurofilament Light Chain) Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purity:	>95% (SDS PAGE)
buffer:	Phosphate Buffered Saline (PBS), pH 7.4
preservative:	0.1% NaN ₃
isotype:	IgG2b
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.