



## PRODUCT SPECIFICATION

catalog number:	<b>A91016</b>
product description:	Calcitonin Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human synthetic Calcitonin
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	BBS 0.1M pH 8.0
preservative:	none
isotype:	IgG1 kappa
concentration:	[lot specific] mg/mL
specificity:	Recognizes natural human Calcitonin
pairing:	This antibody can be used in sandwich ELISA as the conjugate antibody in conjugation with A91018 as the capture
storage:	2-8°C short term
hazard/biohazard:	There is no known biohazard associated with this product.