



## PRODUCT SPECIFICATION

catalog number:	<b>A91066</b>
product description:	Procalcitonin (PCT) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purity:	≥90%
buffer:	10mM PBS, 150mM NaCl
preservative:	0.1% Proclin 300
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
molecular weight:	150kDa
pairing:	This antibody can be used as the capture antibody in combination with A91068 as the detection antibody.
notes/applications:	This antibody's epitope is the KC region and is known to work in CLIA and CMIA.
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.