



PRODUCT SPECIFICATION

catalog number:	A91077
product description:	Procalcitonin (PCT) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant PCT Protein
form:	Liquid
purification method:	Affinity chromatograph purified Fab; 0.2µm Filtered
purity:	>98% by SDS-PAGE
buffer:	PBS, pH 7.2
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL
pairing:	For LF use A91076 as capture with A91077 as conjugate
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.