



PRODUCT SPECIFICATION

CATALOG NUMBER: A91080

PRODUCT:	Anti-PCT (Procalcitonin) Monoclonal Antibody
FORM:	Caprylic acid and ammonium purification from mouse ascites
BUFFER:	10mM sodium phosphate, 0.8% NaCl, pH 7.2 and filtered in 0.2 μ m solution
PRESERVATIVE:	None
STORAGE CONDITIONS:	2-8°C short term; -20°C up to 12 months. Avoid repeated freeze-thaw cycles.
HAZARD/BIOHAZARD:	Standard Laboratory Practices should be followed when handling this material. There is no known hazard or biohazard associated with this product.