



PRODUCT SPECIFICATION

catalog number:	A91093
product description:	Procalcitonin (PCT) Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	≥90% by SDS-PAGE
buffer:	PBS, pH 7.4
preservative:	0.02% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
pairing:	This antibody can be used as the capture antibody in conjunction with A91094 (as the detection antibody).
storage:	2-8°C short term (up to three months); -20°C long term (up to 36 months) avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known hazard/biohazard associated with this product.