



PRODUCT SPECIFICATION

catalog number:	A91140
product description:	Vitamin D Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Vitamin-D derivate-BSA
source:	Murine ascites
form:	Liquid
purity:	> 90%
buffer:	10 mM PBS, 150 mM NaCl, pH 7.4
preservative:	0.1% P300
isotype:	IgG ₁ , kappa
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.