



PRODUCT SPECIFICATION

catalog number:	A91158
product description:	Vitamin D Total (D2 & D3) Monoclonal Antibody
lot number:	[lot specific]
Immunogen:	Vitamin D Derivate - BSA
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1M BBS, pH 8.0
preservative:	0.05% NaN ₃
isotype:	IgG ₁ Kappa
concentration:	[lot specific] mg/mL
Specificity:	Recognizes natural Vitamin D2 and Vitamin D3
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.