



## PRODUCT SPECIFICATION

catalog number:	<b>A91170</b>
product description:	Vitamin D Total (D2 & D3) Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	≥90% by SDS-PAGE
buffer:	PBS, pH 7.4
preservative:	0.02% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
titer:	≥ 1 : 20,000 (ELISA)
storage:	2-8°C short term (up to 3 months); -20°C long term (up to 3 years) avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.