



PRODUCT SPECIFICATION

catalog number:	A91500
product description:	Aggrecan Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Recombinant PG
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	PBS 0.1M pH 7.4
preservative:	0.5g/L NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Recognizes natural human aggrecan. The epitopes recognized are Hyaluronic Acid Binding Region (HARB)
pairing:	This antibody can be used in sandwich ELISA as the conjugate antibody in conjugation with A91501 as the capture
storage:	2-8°C short term
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.