



PRODUCT SPECIFICATION

catalog number:	A91501
product description:	Aggrecan Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Recombinant PG
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	BBS 0.1M pH 8.0
preservative:	none
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Recognizes natural and recombinant human aggrecan.
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjugation with A91500 as the conjugate
storage:	2-8°C short term
hazard/biohazard:	There is no known biohazard associated with this product.