



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A94050**

PRODUCT:	<u>Helicobacter pylori</u> Monoclonal Antibody
IMMUNOGEN:	<u>Helicobacter pylori</u>
SOURCE:	Mouse ascites
PURIFICATION:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
BUFFER:	15 mM Phosphate buffer containing 0.15M NaCl, pH 7.2
PRESERVATIVE:	0.09% sodium azide
SUBCLASS:	IgG <sub>1,k</sub>
STORAGE CONDITIONS:	2-8°C short term; -70°C long term. Avoid freeze-thaw cycles.
APPLICATIONS:	Sandwich ELISA (This antibody can be used as a capture antibody in a <u>Helicobacter pylori</u> ELISA combination with biotinlated, monoclonal antibody A94060489P). Western blot.