



PRODUCT SPECIFICATION

CATALOG NUMBER: A94060

| | |
|---------------------|--|
| PRODUCT: | <u>Helicobacter pylori</u> Monoclonal Antibody |
| IMMUNOGEN: | <u>Helicobacter pylori</u> |
| SOURCE: | Mouse ascites |
| PURIFICATION: | Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography) |
| BUFFER: | 15 mM Phosphate buffer containing 0.15M NaCl, pH 7.2 |
| PRESERVATIVE: | 0.09% sodium azide |
| SUBCLASS: | IgG _{1,k} |
| STORAGE CONDITIONS: | 2-8°C short term; -70°C long term. Avoid freeze-thaw cycles. |
| APPLICATIONS: | Sandwich ELISA (This antibody can be used as a detection antibody in a Helicobacter pylori ELISA combination with the monoclonal capture antibody A94050489P). Western blot. |