



## PRODUCT SPECIFICATION

catalog number:	<b>A94060</b>
product description:	<i>Helicobacter pylori</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	<i>Helicobacter pylori</i>
source:	Murine ascites
form:	Liquid
purification method:	Affinity chromatography
buffer:	15mM phosphate buffer, 0.15M NaCl, pH 7.2
preservative:	0.09% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub> , kappa
concentration:	[lot specific] mg/mL (based on A <sup>1%</sup> <sub>280nm</sub> = 14.0)
pairing:	This antibody can be used in sandwich ELISA as the conjugate antibody in conjunction with A94050 (as the capture antibody).
notes/applications:	Sandwich ELISA
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.