



## PRODUCT SPECIFICATION

catalog number:	<b>A94090</b>
product description:	<i>H. Pylori</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified <i>H. Pylori</i> antigen
source:	Murine ascites
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	PBS, pH 7.0
preservative:	0.02% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A94092 (as the conjugate antibody).
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.