



PRODUCT SPECIFICATION

catalog number:	D-501-C = G83560015
product description:	Sheep IgG (H+L) Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
source:	Donkey
form:	Liquid
purification:	Affinity purified
buffer:	10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% (w/v) NaN ₃
concentration:	[lot specific] mg/mL
specificity:	This antibody reacts with heavy chains on sheep IgG, and with all light chains on all sheep immunoglobulins.
notes/applications:	Serum donors were certified disease free by a registered veterinarian in the U.S.
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.