



PRODUCT SPECIFICATION

catalog number:	D19020								
product description:	T-3 Polyclonal Antibody, Purified								
lot number:	[lot specific]								
immunogen:	T-3 BSA								
source:	Goat								
form:	Liquid								
purification method:	50% (NH ₄) ₂ SO ₄ fraction of whole serum, chromatographically purified (S-100, DEAE) IgG fraction								
buffer:	0.05M PBS, pH 7.4								
preservative:	0.1% NaN ₃								
isotype:	IgG								
concentration:	[lot specific] mg/mL								
affinity constant:	8 x 10 ¹¹								
cross-reactivity:	<table><tr><td>▪ T-3</td><td>100%</td><td>▪ T-4</td><td><1%</td></tr><tr><td>▪ T-2</td><td><1%</td><td>▪ rT3</td><td><1%</td></tr></table>	▪ T-3	100%	▪ T-4	<1%	▪ T-2	<1%	▪ rT3	<1%
▪ T-3	100%	▪ T-4	<1%						
▪ T-2	<1%	▪ rT3	<1%						
storage:	2-8°C								
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.								