



PRODUCT SPECIFICATION

CATALOG NUMBER: D19110

PRODUCT DESCRIPTION:	Rabbit anti-T4
IMMUNOGEN:	T4-BSA
FORM:	Neat Serum
PRESERVATIVE:	None
AFFINITY CONSTANT:	6.0×10^{10}
CROSS REACTIVITY:	T4 100%
STORAGE TEMPERATURE:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.