



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: D24010**

PRODUCT DESCRIPTION:	Rabbit Anti-Testosterone Antiserum	
HOST ANIMAL:	Rabbit	
IMMUNOGEN:	Testosterone-19	
FORM:	Neat serum	
PRESERVATIVES:	None	
CROSS REACTIVITY:	Testosterone	100%
	Danazol	0.004%
	DHT	3.5%
	Nortestosterone	<17.0%
	SHBG	0%
	Other steroids	<0.01%
STORAGE TEMPERATURE:	4-8°C short term, -20°C long term	