



PRODUCT SPECIFICATION

CATALOG NUMBER: D25010

PRODUCT DESCRIPTION:	Rabbit anti Cortisol
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Cortisol-3-BSA
PRESERVATIVES:	None
AFFINITY CONSTANT:	9.0 x 10 ⁹ L/M
TITER:	15,000 test/ml in ¹²⁵ I RIA
STORAGE TEMPERATURE:	2-8°C short term; -20°C long term
CROSS REACTIVITY:	Cortisol 100% 11- deoxycortisol 5.7% Prednisolone 36% Corticosterone 3.3% Cortisone <0.7%
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.