



PRODUCT SPECIFICATION

catalog number:	D41710
product description:	<i>Legionella</i> Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	<i>Legionella pneumophila</i> whole cell prep (ATCC #33152)
source:	Rabbit
form:	Liquid, IgG fraction
purity:	>95%
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
concentration:	[lot specific] mg/mL
titer:	> 1:800 by indirect IFA vs. acetone fixed organisms
specificity:	This polyclonal is known to react with all antigens of the intact microorganism.
cross-reactivity:	This serum has not been absorbed and may react with related microorganisms.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.