



PRODCUT SPECIFICATION

CATALOG NUMBER: D42020

PRODUCT DESCRIPTION:	Anti-Chlamydia trachomatis, IgG Fraction
HOST ANIMAL:	Goat
IMMUNOGEN:	L2+ other serovar groups
BUFFER:	Phosphate Buffered Saline (0.01M, pH 7.2)
PRESERVATIVE:	0.1% sodium azide
TITER:	>1:1,000 by IFA
STORAGE TEMPERATURE:	2-8°C
SPECIFICITY:	Purified elementary bodies, disrupted
UNINFECTED CELL REACTIVITY:	Negative vs. HEp-2 cells and egg yolk sac
CROSS REACTIVITY:	Cross-reacts with Chlamydia psittacii and Chlamydia pneumoniae (TWAR)