



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: D43030**

|                      |  |
|----------------------|--|
| PRODUCT DESCRIPTION: | Rabbit anti- CMV   |
| IMMUNOGEN:           | Purified Cytomegalovirus   |
| FORM:                | IgG fraction from Protein A  |
| PRESERVATIVE:        | None   |
| STORAGE TEMPERATURE: | 2-8°C short term, -20°C long term  |
| SPECIFICITY:         | Reacts against various antigens from all CMV strains. No cross reactivity to other herpes viruses. |