



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D44410**

PRODUCT DESCRIPTION:	Rabbit anti-Hepatitis B Core Antigen (HBcAg)
IMMUNOGEN:	Purified HBcAg
FORM:	BMC-6 carrier polysaccharides, carrier protein
BUFFER:	20 mM tris-borate, 150 mM sodium chloride, pH 7.5
PRESERVATIVE:	0.05% sodium azide
STORAGE TEMPERATURE:	2-8°C, Do not freeze.
REACTIVITY:	This antibody stains predominantly the nuclei of HBV infected liver cells.
STAINING:	This antibody can be used on formalin-fixed, paraffin-embedded tissue sections. Prolonged fixation in buffered formalin can destroy the epitope.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.