



PRODUCT SPECIFICATION

CATALOG NUMBER: D45420

PRODUCT DESCRIPTION:	Rabbit Anti HSV I & II, IgG
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Strain F (human)
BUFFER:	0.01M PBS, pH 7.2
PRESERVATIVES:	0.1% NaN ₃
TOTAL PROTEIN:	4-5 mg/ml
TITER:	>1:2000 by indirect immunofluorescence
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
SPECIFICITY:	ICPS + late structural (virion) antigens
CROSS REACTIVITY:	Cross-reacts with types I & II
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.