



PRODUCT SPECIFICATION

catalog number:	D45500
product description:	HIV p24 Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified native p24 IIIB strain
source:	Goat
form:	Liquid, IgG fraction
purity:	>95%
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
concentration:	4-5 mg/mL
titer:	> 1:1,000 by Western blot
specificity:	This polyclonal sera is specific for p24 (as shown by predominant band in Western)
cross-reactivity:	This serum does not cross-react with human B or T cells and is not reactive with human or bovine sera.
storage:	2-8°C (<6 months); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.