



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D46000**

PRODUCT DESCRIPTION:	Goat anti RSV, IgG Fraction
HOST ANIMAL:	Goat
IMMUNOGEN:	Human RSV isolate, confirmed
PURITY:	>95%
BUFFER:	Phosphate buffered saline (0.01M, pH 7.2)
PRESERVATIVE:	0.1% sodium azide
TITER:	>1:480 by IFA
STORAGE TEMPERATURE:	2-8°C
SPECIFICITY:	All viral antigens
UNINFECTED CELL REACTIVITY:	Negative vs. HEp-2 and WI-38 cells
CROSS REACTIVITY:	Does not cross-react with Para 1-3, Influenza A & B or Adenovirus by IFA
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.