



PRODUCT SPECIFICATION

CATALOG NUMBER: D46530

PRODUCT DESCRIPTION:	Pseudomonas sp. Polyclonal
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Whole cells of Pseudomonas aeruginosa
FORM:	Neat antisera
PURITY:	>95%
PURIFICATION:	Purified IgG fraction
BUFFER:	0.01M phosphate saline buffer, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
SPECIFICITY:	Soluble + structural antigens of whole bacterium.
CROSS REACTIVITY:	Antiserum is not absorbed and may crossreact with related microorganisms.