



PRODUCT SPECIFICATION

catalog number:	D48030
product description:	<i>Streptococcus pneumoniae</i> Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Whole cell blend of multiple serotypes
source:	Rabbit
form:	Liquid, IgG Fraction
purity:	>95%
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
concentration:	[lot specific] mg/mL
specificity:	This polyclonal is known to react with types 3, 4, 6, 7, 9, 14, 18, 19 and 23.
cross-reactivity:	This serum has not been absorbed and may react with related microorganisms.
storage:	2-8°C (six months or less); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.